

**Hiflex<sup>5</sup> the perfect complement for your therapy furniture – simply modern and functional**



# Hiflex<sup>5</sup> – wherever it is needed



The Hiflex<sup>5</sup> can be moved easily in every desired position and the base frame with its small shell matches the Bionic Therapy Chairs perfectly.



The multifunctional table can be height adjusted smoothly and fast; the big tray is tiltable in two steps.



There is a practical hook made of Aluminium for clothes or bags on the inside of the column. The Hiflex<sup>5</sup> can be equipped with a useful container or...



..., if desired, with a big newspaper holder. Two separate trays offer much space for medical equipment or personal belongings during treatment.

With the simple and continuous height adjustment by foot pedal the Hiflex<sup>5</sup> meets your patients' wishes. Additionally, the table top is tiltable up to 30° in two steps.

Four smooth-running wheels as well as a practical handle allow an easy movement and positioning of the Hiflex<sup>5</sup>.

At the inner side of the column there is a stable hook made of Aluminium for storage of bags or clothes. Alternatively, the Hiflex<sup>5</sup> is equipped with a container or a newspaper holder and therefore offers enough space.

The Hiflex<sup>5</sup> is available in light grey (RAL 7035) to match the Bionic Therapy Chairs perfectly.

## Technical Data:

Dimensions	
Total width	40 cm
Total length	78 cm
Height adjustment	74 – 97 cm
Removable Trays	
Big tray	55 x 39 cm
Small tray	15 x 39 cm
Safety	
Safe maximum load	15 kg
Optional Accessories	
Container, H x L x W	22 x 20 x 30 cm
Newspaper holder, H x L x W	22 x 13 x 31 cm
Available colours	
Light grey	RAL 7035



Bionic Medizintechnik GmbH ... a JMS-corporation · Max-Planck-Straße 21 · D-61381 Friedrichsdorf/Germany  
Phone: +49(0)6172-75760 · Fax: +49(0)6172-757610 · info@bionic-jms.com · www.bionic-jms.com



**TekMed™**

1300 720 727



12 Bellevue Crescent Preston Victoria 3072 Australia  
T (03) 9487 9999 F (03) 9484 8792  
info@tekmed.com.au www.tekmed.com.au

**TekMed NZ™**

Ph: 0800 456 217

info@tekmed.co.nz

www.tekmed.co.nz